



NEW ACCOUNT SETUP:

1. Contact representative below or visit our website for registration
2. Complete the application (with qualified NPI number)
3. Submit application with customer service or upload online

ORDER PROCESS:

Contact representative below or visit our website for more information.

Please be aware of the following:

- Overnight deliveries contain dry ice
- Delivery available Tuesday - Saturday
- Next day shipping is available for orders placed by 2:00PM Central Time

Distributed by:



MISSION & VISION

Signature Biologics uses its innovative techniques to manufacture human placental-derived products to support and improve the natural healing processes of the body. We strive to produce best-in-class products that better the quality of life of patients.

Our vision is to become the premier, trusted source of innovative placental-derived therapeutic solutions in the nation, characterized by the highest quality products and improved functional patient outcomes.

VALUES

QUALITY: Compliance with the highest regulatory standards to ensure products consistently meet established specifications.

INTEGRITY: To conduct all business operations with the utmost commitment to transparency, trustworthiness, and golden rule ethics.

INNOVATION: To utilize a constant focus on research to develop novel and unique processes and best-in-class products.



11496 Luna Rd, STE 800
Dallas, TX 75234
972-807-5916
info@signaturebiologics.com
www.signaturebiologics.com

PR-SC 001 rv.04 10MAR2020



Signature CORD

Human Umbilical Cord Tissue

www.signaturebiologics.com



Intended Use:
Cushioning Agent

Available Sizes in
Vial Volume (Tissue Weight):
1.0mL (75mg) & 2.0mL (150mg)



TISSUE PRODUCT

DESCRIPTION:

Cushioning agent derived from Wharton's Jelly tissue for homologous supplementation.

Our cord product has been processed through a 500 (Standard) or 200 (Micro) micron filter.


MEASUREMENT:

Signature Cord 2x = 75mg tissue per mL

DELIVERY:

Signature Cord can pass through a 21-gauge or larger diameter needle.

Signature Cord Micro can pass through a 25-gauge or larger diameter needle.



*Actual image of Signature Cord tissue product.
Color added for marketing purposes only.

USE & PROPERTIES



NATURAL BIOLOGIC ALTERNATIVE

The umbilical cord represents the link between mother and fetus. Our Signature Cord product consists of isolated Wharton's Jelly tissue and does not contain HLA antigens that are normally associated with the immune response.



STRUCTURAL TISSUE PRODUCT

Using our proprietary Signature Processing techniques, our tissue maintains the original relevant characteristics relating to its utility for soft tissue supplementation.



MINIMALLY INTRUSIVE

Less invasive and less painful than adipose tissue extraction and bone marrow aspiration-based therapies.



CRYOPRESERVED AS A FROZEN LIQUID IN AN EASY-TO-OPEN VIAL

Cryopreserved in 10% DMSO for storage below -65 °C. Thaws in gloved hand in less than 5 minutes for application.

MARKET & COMPLIANCE



OUR PRODUCTS ARE DERIVED ONLY FROM HEALTHY, PRE-SCREENED LIVE C-SECTION BIRTHS.



TO ENSURE PATIENT SAFETY, OUR DONOR SELECTION, UMBILICAL CORD RECOVERY AND PROCESSING PROCEDURES MEET OR EXCEED ALL APPLICABLE INDUSTRY STANDARDS.



ALL PRODUCTS UNDERGO INDEPENDENT STERILITY TESTING AND ARE MANUFACTURED IN A GTP FACILITY.



SIGNATURE CORD IS (HCT/P) COMPLIANT UNDER SECTION 361 OF THE PUBLIC HEALTH SERVICE ACT ACCORDING TO 21 CFR PART 1271.10. IT IS REGULATED AS A HUMAN CELL AND TISSUE PRODUCT - FOR HOMOLOGOUS USE ONLY.



PRESCRIPTION USE ONLY FOR SALE TO LICENSED HEALTH CARE PROFESSIONALS.

